











Ocular Vitrectomy Lens Rings

	<p>Journal Reference: Ophthalmic Surgery and Lasers Vol. 27, No. 10, p. 891, October 1996</p>	<p><i>Landers Rings Designed with: Maurice B. Landers, III, M.D., Chapel Hill, NC</i></p> <p><i>Tano Ring Designed with: Yasuo Tano, M.D., Osaka, Japan</i></p> <p><i>Foxman Ring Designed with: Brett Foxman, M.D., Northfield, NJ</i></p>
---	---	---

Ocular Vitrectomy Lens Rings are designed to be used with Ocular vitrectomy lenses in both the Landers Lens Ring System (OLVS-3 or -3N) and the HRI Lens Set (OLVS-HRI). OPV-R is designed for use with the pediatric sized vitrectomy lenses in the Ocular Pediatric Vitrectomy Lens Set (OPV-S).

	Product Code/ Lens Name	Design
	<p>OFV-4 Foxman Vitrectomy Lens Ring</p>	<p>4 down struts. Requires no sutures. Designed for use with a 25g cannula</p> <p>Technique</p> <p>§ Place ring on eye.</p> <p>§ Mark sclera for the two superior cannulas using the 3 & 4mm markings on the ring.</p> <p>§ Remove ring. Insert cannulas.</p> <p>§ Replace ring.</p>
	<p>OLV-1 Landers Vitrectomy Lens Ring <i>(Included in set OLVS-3)</i></p>	<p>Stainless steel ring with two suture down struts.</p>
	<p>OLV-1-TN Landers Tall Notched Vitrectomy Lens Ring <i>(Included in set OLVS-3N)</i></p>	<p>This stainless steel ring is centered on the cornea. Three notches are designed in the top of the ring for suture placement in the sclera. Height is 3.2mm.</p>
	<p>OLV-1-IR Landers Irrigating Vitrectomy Lens Ring</p>	<p>Stainless steel ring features an irrigation port. Sutures secure the two struts to the sclera, which allows blood to be irrigated away and keeps the cornea moist.</p>
	<p>OLV-1-IN Landers Notched Irrigating Lens Ring</p>	<p>Stainless steel ring features an irrigation port. Three notches for sutures in the top of the ring. Height is 3.2mm.</p>
	<p>OLV-1S Landers Silicone Ring</p>	<p>Can be used with any Landers System Lens. It allows the surgeon to change lens positions to obtain the optimum viewing angle. It is useful when using a prismatic lens for peripheral vitrectomy procedures. (4 per pack)</p>
	<p>OLV-1-4P Landers 4 Post Vitrectomy Lens Ring</p>	<p>The Landers 4 Post Vitrectomy Ring consists of a stainless steel ring with an inside diameter of 11.5mm and four posts set at a 30° angle. The posts are 2mm tall and the ring is 3.2mm tall. The new design permits easier scleral depression in the 9:00 and 3:00 regions. Two posts, with space in between, allow compensation for slightly misplaced sutures to easily maintain ring centration on the cornea.</p> <p>Technique</p> <p>§ The Landers Ring is positioned with the cornea geometrically centered within the ring.</p> <p>§ Using a 5.0 Mersilene suture, the suture needle is inserted immediately next to the limbus at the 3:00 and 9:00 position using a single armed suture.</p> <p>§ The suture is then cut midway between the 3:00 and 9:00 posts.</p> <p>§ The suture is then tied around both posts on each side.</p> <p>§ If the sutures are not placed precisely at the 3:00 and 9:00 positions, the ring may become decentered, compromising the view of the retina.</p> <p>§ The first suture is tied loosely holding the ring in place. <i>If tied too tightly, the ring will tilt upward and may even flip over.</i></p> <p>§ The second suture is tied tightly in place.</p> <p>§ Finally, the first suture is re-tied tightly in place.</p> <p>§ If the ring is slightly off center, the knot can be tied over one of the posts rather than both, thus maintaining the centration.</p>

	OPV-R Pediatric Vitrectomy Lens Ring	Smaller version of the OLV-1. <i>(Included in Pediatric Set OPV-S)</i>
	OTN-R Tano Vitrectomy Lens Ring	4 down tabs. Requires only 1 circumferential suture allowing for easy repositioning and removal.

Cleaning & Disinfection

See Cleaning Method 3



2255 116th Ave NE, Bellevue, Washington 98004-3039 USA
 T: 425-455-5200 or 800-888-6616 F: 425-462-6669
 E: ocular@ocularinc.com I: www.ocularinc.com

© 2001 Ocular Instruments
 597113275