


Ocular Indirect Lens Sterilizing Tray



PRODUCT
CODE

OI-ST

CE

Design

§ Sterilizing tray for MaxLight[®], MaxField[®], and MaxAC[®] lines of Ocular Indirect Ophthalmoscopy lenses.

§ It makes ethylene oxide and cold sterilization of lenses quick and easy.

§ Use for steam autoclaving of MaxAC[™] lenses.

§ This stackable, durable tray gently holds lenses to protect them during sterilization.

Preparation of Lens(es)

§ Clean with alcohol wipe.

§ Place in Indirect Lens Sterilizing Tray.

§ Proceed with sterilization instructions.

